

Test 38: Warning lights On acceptance & Daily


Radiation warning light at entrance of all x-ray rooms, excluding theatres

Must work when beam is activated

| TEST NO.38: RADIATION WARNING LIGHT - RM _____ | | | | |
|--|------|------|------|-------------------|
| DATE | NAME | PASS | FAIL | CORRECTIVE ACTION |
| | | | | |
| | | | | |
| | | | | |



Warning lights



Warning lights for General, CT & Fluoroscopy units (excluding theatres):

A **red** warning light, which is only activated when the beam is on and when fluoroscopy is in progress, must be mounted in a conspicuous place outside the entrance to the x-ray rooms

- Required warning lights must be in working order

