

## Digital Mammography

### Test 56-66: Repeat and reject analysis 3 monthly

For this section use European guidelines for quality assurance in breast cancer screening and diagnosis (Ref 4)

Repeat and reject analysis

**AEC device** On acceptance & weekly

Long term reproducibility

**Image receptor homogeneity** On acceptance & Weekly

## Digital Mammography

### Test 67-68: Small field Digital Mammography

### Tests 69-74: Mobile mammography systems

Must ensure that all freely moveable objects/equipment are **firmly locked** or strapped down before moving.

Perform visual check of breast support and associated equipment for possible damage.

#### After moving check:

- Compression device

- Alignment of x-ray beam to image receptor

- AEC system

- Image quality