

Test 3: Appropriate technique chart displayed at x-ray unit

6 monthly

SITE: _____				
TEST NO. 3: VISUAL CHECK OF TECHNIQUE CHART - ROOM__				
DATE	NAME	PASSED	FAILED	CORRECTIVE ACTION

Do not copy charts from other rooms.
Check synergy with EI/LGM/S value.

ONE size exposure chart does **NOT** fit all!

BELIGTINGSKAART / EXPOSURE CHART

ONDERSOEK EXAMINATION	POSISIE POSITION	CM PA/AP	CM LAT	FOKUS FOCUS	FFD	ROOSTER BUCKY	kV	mAs
CHEST								
female	PA	18-20		BF	180	JA	109	2.5
female	PA	20-23		BF	180	JA	109	2.8
female	PA	24-26		BF	180	JA	109	3.2
female	LAT		25-27	BF	180	JA	117	4.0